





Scanned + Archived
By:
FindMyManual.org



EYE SPA

USER MANUAL



PatriotHealthAlliance

2920 Berry Hill Dr., Suite 300
Nashville, TN 37204

(800) 230-8956 | PatriotHealthAlliance.com

V1.0

CONGRATULATIONS! Your Eye Spa is the perfect antidote to a stressful day, and will help you relax from all the pressures and busyness of everyday life.

The Eye Spa is an electric eye/temple area massager that uses vibrational, layered compression with heat to ease tension and deliver relaxation to tired, stressed muscles.

Use your Eye Spa at home, work, or take it with you on your travels for a rejuvenating experience almost anywhere.

INSTRUCTIONS FOR USE

1. Take your Eye Spa out of its packaging.
2. Find a relaxing spot to use your Eye Spa, such as a couch or bed.
3. Locate the on/off button. Press longer to turn on, and shorter to switch between modes: dynamic, relaxation, soothing, or sleep. Choose between 10-20 minute cycles.
4. Optional Bluetooth connection can play music from your phone or tablet.
5. The internal battery is rechargeable with the included USB cable. When running low, the battery charge indicator light will flash. Simply recharge as needed to continue enjoying your Eye Spa.

FEATURES AND SPECS

- USB-C Input: 5V/1A
- UV Output: 2W
- Wireless Output: 5W
- Color: White
- Materials: ABS Plastic
- Dimensions: 8" x 4.3" x 2.7"
- Weight: 0.71lbs

Warning: The Eye Spa is not recommended for individuals with any kind of eye disease, or recovering from surgery in the eye or head area. Persons who are at risk of, or who have suffered a stroke, and pregnant women are also strongly discouraged from using Eye Spa.

Eye Spa is not a toy. Do not give this device to children.