BOLT +

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INTRODUCTION


MUSCLE RECOVERY WITH JOLT™ BOLT+

The JOLT™ Bolt+ Muscle Massage Gun is a portable, rechargeable vibration massage device. JOLT™ Bolt+ helps to increase blood flow and effectively alleviates the discomfort caused by the accumulation of lactic acid after sports and thus reduce discomfort for people suffering from muscle stiffness / aches.

JOLT™ can be used in the following cases:

- For warming up before a physical effort
- For the recovery after a sporting effort
- To relieve aches and joint stiffness
- To stimulate blood circulation
- To relieve stress and tension
- Massage of knots and trigger points
- Suitable for beginners, amateurs and professionals
- To relax and to loosen the tissues
- To reduce lactic acid build up in the muscles
- To increase joint mobility
- To fluidize edema in the swollen area
- To relieve aches and pains
- To relieve tendinitis


CONTENTS

- 1 1 JOLT Bolt+ massage gun
- 2 6 massage heads
- 3 1 DC charger cable
- 4 2 manuals (Massage program and technical instructions)
- 5 1 carry case

GETTING STARTED GUIDE

1 Description of the JOLT™ BOLT+

2 ***JOLT™ BOLT+ steps for use***

Our massage gun has **5 classic speeds** and **3 intelligent speeds**:

5 classic speeds

(lit in **blue** above the control button): 1500rpm / 1900rpm / 2300rpm / 2700rpm / 3200rpm (rpm = repetitions per minute)

3 Smart-Hit speeds

(illuminated in **orange** below the control button):

Increasing speed (muscle activation before exercise)

Alternating speed (relaxation scenarios)

Decreasing speed (muscle relaxation after exercise)

1. Turn on the unit using the ON/OFF button located under the handle.

2. Press the center button on the control interface to start the gun in conventional speed 1 (**blue color**).

Press a second time to switch to speed 2 and so on up to speed 5.

3. To switch to the smart speeds, press the center button for 1-2 seconds. There are 3 different smart speeds and they are illuminated in **orange**.

4. To stop and turn off the gun:

- If you are using the «classic» speed mode: press the central button on the control interface repeatedly until you reach speed 5, then press it one last time. Next, switch the button under the handle to the OFF position.
- If you are in «smart» speed mode, press the center button on the control interface until you reach smart speed 3, then press it one last time. Then turn the button under the handle to the OFF position.
- Or turn the button under the handle directly to the OFF position without following the previous steps.

5. The battery charge status is displayed at the bottom of the handle. When you want to put the massage gun on charge, turn it off and plug it in, you will see the % displayed and gradually increase until you reach 100% (full and finished charge).
6. The unit has a timing feature and will automatically shut off after 15 minutes of continuous operation (for safety to avoid overheating and to preserve the integrity of your muscles). You can turn it back on after 1 minute of pause.
7. When the massage head is pressed too hard, the unit will start its automatic protection function and stop. It can be used normally after the unit is turned on again.

3 *Description of the JOLT™ BOLT+ Massage Heads*

With the massage gun turned off, slowly insert the massage head into the barrel until it stops. As shown in the illustration, the position of the red line on the massage head should be aligned with the position of the red line on the gun. To remove the massage head, you just need to pull on the head to remove it.

JOLT™ Bolt+ can be used with 4 massage heads depending on the parts of your body you want to massage.

Round head

Its spherical design with soft coating allows to work large muscle groups such as the back of the thigh and the calf, the triceps, ... It is mainly suitable for new users who are sensitive to vibrations.

Flat head

It allows you to relax any muscle of the body for a classic massage.

Cylindrical head

The conical shape of this head allows an impacting and toning massage with a significant pressure depth. It is suitable for the meridians, soles of the feet, palms of the hands, Achilles heel, rhomboid muscles of the back, ...

U-shaped head

Its shape fits perfectly the spine and cervical/neck. It can also fit the shoulder and other bilateral muscles (do not touch directly the spine and cervical to avoid injury).

Round head XL

Big brother of the «classic» round head, it has the same characteristics, but its larger size allows you to train larger muscles: Hamstrings, quadriceps, glutes.

Cushioned head

This head is a variation of the flat head, we have added on its top a soft layer to adapt to the roughness of your body. It is suitable for soft tissues (tissues that support the structures and organs of the body such as muscles), internal and external muscles of the abdomen, chest and legs. The air it contains cushions the percussions for absolute comfort.

PRECAUTIONS

1 *Warning*

- Do not place the device in water, keep it away from liquids or heat sources, and prevent water from entering the product's slots. They must be free of dust and debris . Store it in a dry place.
- Do not drop or use the product incorrectly.
- Do not remove the screws or attempt to disassemble it.
- Do not change or modify the product in any way.
- Do not use the device continuously for more than 1 hour and let it stand for 30 minutes before using it again.
- Only charge the device with the charging cable. Charge the device every 6 months if you are not using it for a long time
- Unplug the charger from the device after charging or before use.
- Do not use or charge the device unattended.
- Use a soft cloth for daily maintenance. Do not use engine oil, organic solvents or others.

2 *Health warning*

- Do not dispose of this product in a conventional trash can at the end of its life. Dispose of it in a recycling garbage can or return it to the JOLT™ Service Center. We will recycle it in an environmentally friendly way.

- The use of JOLT™ is strongly discouraged in the following cases:
 - Pregnancy
 - Diabetes with complications such as neuropathy
 - Retinal damage
 - Heart disease, with pacemaker or defibrillator
 - Recent surgery (less than 90 days)
 - Epilepsy
 - Migraines
 - Joint prostheses
 - Microtrauma
 - Stenosis of the blood
 - Tumor area
 - Sensitive areas: head, face, cervical vertebrae, spine
 - Positions close to the bone with less muscle coverage: humerus, back of the foot, back of the hand
 - Other health concerns you may have.
- To reduce the risk of personal injury, this product should be used as follows:
 1. Use is exclusively for adults, use only on dry parts of the body on clean skin.
 2. Use the product only on soft tissues of the body (without pain or discomfort), do not use on hard parts like skull or body bones.
 3. Only use massage heads which are recommended for the muscle groups involved.
 4. Whether controlled by deep or surface pressure, too much massage can cause hematoma; check the massage area regularly and stop using it as soon.

3 Environmental protection

Do not dispose of this product in a conventional trash can at the end of its life. Dispose of it in a recycling garbage can or return it to the JOLT™ Service Center. We will recycle it in an environmentally friendly way.

TECHNICAL SPECIFICATIONS

- Model : JOLT BOLT+
- 5 classic speeds from 1500 to 3200rpm (repetitions per minute)
- 3 intelligent speeds: increasing speed / decreasing speed / alternative speed
- Digital display
- 6 interchangeable massage tips (included)
- Brushless motor 55m & 125W
- Average battery life of 5-7h depending on the selected power level
- Charging time : 3h
- Rechargeable Lithium battery 2500mAh
- Charging with DC charger (included)
- Pressure from 20Kg
- Quiet Technology lower than <50 decibels
- Amplitude : 14mm
- Comfortable non-slip and anti-sweat grip
- On/Off indicator
- Speed indicator
- Charging indicator
- 990g net weight
- Massage gun size: 25cm x 18cm x 6cm
- Storage box included

FAQ & WARRANTY

Any questions or particular concerns?

Please visit our Frequently Asked Questions (FAQ) page for answers to all your questions:

*If you cannot find a solution, please contact our customer service at the following address: **hello@getjolt.fr***

MORE QUESTIONS / ANSWERS

The warranty period for this product is 24 months from the date of purchase. Wear and tear due to appearance and mechanical operation are not covered by the warranty.

The warranty does not cover:

1. Accidental damage, misuse, or transportation.
2. Improper maintenance.
3. Operation without following the manual instructions.
4. Unable to provide a formal purchase certificate.

JOLT

www.getjolt.fr