Össur Americas

27051 Towne Centre Drive
Foothill Ranch, CA 92610, USA
Tel: +1 (949) 382 3883
Tel: +1 800 233 6263
ossurusa@ossur.com

Össur Canada

2150 – 6900 Graybar Road
Richmond, BC
V6W 0A5, Canada
Tel: +1 604 241 8152

Össur Europe BV

De Schakel 70
5651 GH Eindhoven
The Netherlands
Tel: +800 3539 3668
Tel: +31 499 462840
info-europe@ossur.com

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
ÖSSUR®

LIFE WITHOUT LIMITATIONS

Scanned + Archived
By:
FindMyManual.org

Instructions for Use


UNLOADER ONE® X

Össur hf.

Grjótháls 1-5
110 Reykjavík
Iceland

		4
EN	Instructions for Use	6
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	18
NO	Bruksanvisning	21
DA	Brugsanvisning	24
SV	Bruksanvisning	27
FI	Käyttöohjeet	30
NL	Gebruiksaanwijzing	33
PT	Instruções de Utilização	36
PL	Instrukcja użytkowania	39
CS	Návod k použití	42
JA	取扱説明書	45
ZH	中文说明书	48
KO	사용 설명서	51

• **Figure 3:** Ask user to sit down and extend leg while positioning the brace. The Quick Fit buckles and corresponding keyholes are color matched. Place lower buckle button into calf shell keyhole (blue) and use the palm of the hand to snap the buckle closed. Position calf strap above widest part of gastroc, fold strap at appropriate length into alligator clip and fasten strap so that it keeps brace securely and correctly positioned on leg.

• **Figure 4:** Place upper buckle button into thigh shell keyhole (yellow) and use the palm of the hand to snap the buckle closed. Adjust thigh strap to appropriate length by tensioning around leg and folding into alligator clip.

• **Figure 5:** With both SmartDosing dials at lowest setting, adjust the upper DFS strap length until it sits firmly against the leg, and then fold it into alligator clip. At this point, user should not experience any tension or unloading. Adjust the lower DFS strap in the same way.

• **Figure 6:** Ask user to bend knee with the foot flat on the floor.

• **Figure 7:** Turn upper and lower SmartDosing dials clockwise until indicators are in the middle slot.

- Have user stand up and take a few steps to verify correct positioning of the brace and tightness of the straps.
- Determine optimal DFS strap tension based on user's pain relief feedback.
- If user requires more or less tension with indicator at mid setting, adjust the length of the DFS straps accordingly.
- Aim for final SmartDosing setting at middle indicator slot since this will give user ability to adjust dosing during daily life activities.

• **Figure 8:** When final fit has been confirmed, trim the straps to the appropriate length starting with the calf strap so that brace sits correctly on leg while trimming other straps.

• **Figure 9:** With user in standing position, insert fingers into popliteal pad openings and move to an ideal position in popliteal area.

BRACE REMOVAL

1. Ask user to sit with knee flexed and foot flat on floor.
2. Open both upper and lower buckles.
3. Pull buckle buttons out of keyholes.










CLEANING

Insert entire brace into designated washbag, without removing pads or straps.

Machine wash using a gentle wash cycle.

Wash with a gentle detergent.

NOTES:

- Avoid contact with chlorinated water or salt water. In case of contact, rinse well with tap water.
- Do not wash using bleach as this will affect durability and may cause skin irritation.
- Do not wash using fabric softener.
- Do not tumble dry.

PRODUCT MATERIAL LIST

Aluminum, silicone, nylon, stainless steel, acetal, polypropylene, thermoplastic elastomer, polycarbonate, polyester, polyurethane, spandex, lithium grease, acrylic adhesive.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

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