
Scanned + Archived
By:
FindMyManual.org

medela 

Swing™

EN Instructions for use | ES Instrucciones de uso | FR Mode d'emploi

medela 

www.medela.com


Medela, Inc.,
1101 Corporate Drive
McHenry, IL 60050, USA
Ph/Tel./Tél.: (800) 435-8316 / (815) 363-1166
Fax/Téléc.: (815) 363-1246
Email/Correo electrónico/Courriel:
customer.service@medela.com

Medela Canada, Inc.,
4160 Sladeview Crescent
Unit #8, Mississauga, Ontario,
Canada L5L 0A1
Ph/Tel./Tél.: (800) 435-8316
Fax/Téléc.: (800) 995-7867
Email/Correo electrónico/Courriel:
info@medela.ca

Printed in the USA. Medela and 2-Phase Expression are registered trademarks of Medela Holding AG.
Swing, Quick Clean, Micro-Steam, Maximum Comfort Vacuum and PersonalFit are trademarks of Medela, Inc.

Impreso en los Estados Unidos de América. Medela y 2-Phase Expression son marcas registradas de Medela Holding AG.
Swing, Quick Clean, Micro-Steam, Maximum Comfort Vacuum y PersonalFit son marcas comerciales de Medela, Inc.

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Important Safeguards

When using electrical products, especially when children are present, basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

DANGER: To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place or drop into water or other liquid.
- Do not reach for a product if it has fallen into water. Immediately unplug from electrical outlet.

WARNING: To avoid fire, electrocution, or serious burns:

- Do not leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when this product is used near children or persons with disabilities.
- Use the product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep cord away from heated surface.
- Never use while sleeping or when drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors, or operate where aerosol (spray) products are being used where oxygen is being administered.
- Always inspect power adaptor and battery pack wires prior to use for damage or exposed wire. If damage is found, immediately discontinue use of power adaptor or battery pack and call Medela Customer Service at 1-800-435-8316.

The warning symbol identifies all instructions that are important to safety. Failure to observe these instructions can lead to injury or damage to the breastpump. When used in conjunction with the following words, the warning symbols stand for:

WARNING Can lead to serious injury or death.

CAUTION Can lead to minor injury.

NOTE Can lead to material damage.

INFO Useful or important information that is not related to safety

WARNING: To avoid health risk and reduce the risk of injury:

- This product is intended for use by a single user only and should not be shared between users.
- Do not drive while pumping.
- Never use while pregnant, as pumping can induce labor.
- Clean and sanitize all parts that come in contact with your breast and breastmilk prior to first use.
- Wash all parts that come in contact with your breast and breastmilk after every use.
- Inspect all appropriate pump components before each use.
- Do NOT continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- Do not thaw frozen breastmilk in a microwave or in a pan of boiling water.
- If tube becomes moldy, discontinue use and replace tube.
- If breastmilk backup occurs you must clean tube before your next pumping session.

CAUTION: Can lead to minor injury:

- Do NOT wrap cord around adaptor body.
- Use only the power adaptor that comes with the Swing Breastpump.
- Use only the Battery Pack that comes with the Swing Breastpump.
- Pump only with the breastpump in an upright position.
- Do NOT shorten tubing.
- Make sure the voltage of the power adaptor is compatible with the power source.
- Plug the power adaptor into the breastpump first and then into the wall socket.
- Do NOT use antibacterial or abrasive cleaners/detergents when cleaning breastpump or breastpump parts.
- Never put breastpump motor unit in water or a sterilizer, as you can cause permanent damage to the breastpump.
- Do not attempt to remove the breastshield from your breast while pumping. Turn the breastpump off and break the seal between your breast and breastshield with your finger, then remove breastshield from your breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breastshield with your finger and remove the breastshield from your breast.
- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breastmilk if bottles or components become damaged.

SAVE THESE INSTRUCTIONS

Questions? Visit www.medela.com or call us at 1-800-435-8316.

Table of Contents

1. Cleaning	8-13
1.1 Before using for the first time	8
1.2 After each use	8-9
1.3 Sanitize daily	10-11
1.4 Tubing and motor unit cleaning	12-13
2. Assembly.....	14-21
2.1 Assembly of pump kit	14-17
2.2 Correctly fitting breastshields	18
2.3 Powering your pump	19-21
2.3.1 Power adaptor care	19
2.3.2 Battery compartment.....	20-21
2.3.3 Traveling outside of the U.S. or Canada	21
3. Pumping.....	22-26
3.1 Device operation.....	22-25
3.2 Mobility.....	26
4. Storing, Preparing & Feeding Breastmilk	27
4.1 Storing breastmilk.....	27
4.2 Preparing breastmilk.....	27
4.3 Feeding breastmilk.....	27
5. Troubleshooting	28-29
6. Supplemental Information/Disposal Information.....	30
7. Warranty.....	31

Product Description

Take a moment to read through this entire instruction manual before using this product for the first time.

PLEASE SAVE THESE INSTRUCTIONS.

It is best to wait until your breastfeeding routine is established (approximately 4 weeks) before expressing breastmilk, unless otherwise advised by your healthcare professional.

Glossary of Terms

2-Phase Expression® technology – research based technology that mimics a baby's natural nursing rhythm.

Expression Phase – slower sucking/pumping rhythm for gentle and efficient milk removal as quickly as possible.

Maximum Comfort Vacuum™ – highest vacuum setting where a mother still feels comfortable during pumping. It is different for every mother.

Stimulation Phase – fast sucking/pumping rhythm to stimulate the milk ejection reflex and to start the milk flowing.

Intended Use

For lactating women to express and collect milk from their breasts to complement breastfeeding.

If you have medical reasons or other needs for exclusively pumping, it is recommended that you use a hospital grade breastpump such as our Symphony breastpump. To learn more, visit us at www.medelabreastfeedingus.com.

Product Description

Swing™ is a personal use electric breastpump that includes 2-Phase Expression technology and is capable of single pumping.

Product Description

1x PersonalFit™

Breastshield M (24mm)
(for other breastshield sizes
see **Section 2.2**, page 18)
87073 (24 mm)

1x Pump Body

8007196

2x Valves

87089

2x Membranes

87088

1x Bottle

6007135

1x Bottle Stand

8100462

1x Neck

Strap

8000734

1x Swing Storage Bag

8000838

1x Swing Tubing

8007215

Overview of LEDs

lights up	On
does not light up	Off
even, regular flashing	Stimulation phase
light constantly on	Expression phase
rapid flashing	Min./Max. vacuum

LED

On/Off button

Reduce
vacuum

Increase
vacuum

Let-down
button

1x Motor Unit

0300011

Belt clip

Tubing port

Battery
compartment

Power
adaptor port

1x Backplate

2x Lids

8107183

1x Power Adaptor


9200043

1. Cleaning

1.1 Before using for the first time

Clean and sanitize the parts below before using your pump for the first time. Follow the cleaning and sanitizing instructions in **Sections 1.2 (start at step 4) and 1.3** that follow.


Caution


- Only use drinking-quality tap or bottled water for cleaning.
- Take apart and wash all parts that come in contact with the breast and breastmilk immediately after use to avoid dry up of milk residues and to prevent growth of bacteria.

1.2 After each use


Clean all parts that come in contact with your breast and breastmilk. See image in **Section 1.1**.

Separate breastshield from pump body.

Twist and pull the yellow valve off the pump body.

Remove the white membrane from the yellow valve.

Washing Options:

Washing in sink:


- Soak all separated parts in warm soapy water for 5 minutes.
- Clean with a clean dish-cloth or soft brush.
- Rinse all separated parts with clear water.
- Allow all parts to air dry in a clean area.
- Store dry parts when not in use.

or

Washing in dishwasher:

- Wash all separated parts on top rack of dishwasher.
- Allow all breastpump parts to air dry in a clean area.
- Store dry parts when not in use.

Questions? Visit www.medela.com or call us at 1-800-435-8316.

1. Cleaning (cont.)


! Note

- **Wash** breastpump kit parts after each use.
- **Sanitize** breastpump kit parts once daily.
- The pump kit requires no maintenance.
- Take care not to damage parts of the pump kit during cleaning.
- If the individual components of the pump kit are cleaned in a dishwasher parts may be discolored by food pigments. This has no impact on the function.
- Store the pump kit in a clean bag/container until next use.

1.3 Sanitize daily


Wash hands thoroughly.


Refer to image in Section 1.1


Separate all parts that come in contact with your breast and breastmilk.

- Fill a pot with enough water to cover all parts.
- Bring water to a boil.
- Place parts in boiling water for 10 minutes.

Allow water to cool and gently remove parts from water with tongs.

or

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts when not in use. Do NOT store wet or damp parts.

You may sanitize your kit parts by following the instructions on Medela's Quick Clean™ Micro-Steam™ bag.

1. Cleaning (cont.)

Caution

- If tubing becomes moldy, discontinue use and replace tubing. Contact Medela Customer Service at customer.service@medela.com or 1-800-435-8316.
- Do not wash power adaptor.

Note

- Tubing should be washed if dirty or milk is present. Follow instructions in **Section 1.4**.


1.4 Tubing and Motor Unit Cleaning

1

Turn off breastpump.

2


Unplug breastpump from power source.

3


Remove tube from motor unit.

4

Remove tube from pump body.

5

Rinse tube in cool water to remove breastmilk.

6

- Wash tube in warm soapy water.
- Rinse tubing with clear water.
- Hang to air dry.


7

Wipe off with a clean, damp cloth.

2. Assembly

Parts you will need for this step:


! Note

- Inspect valves for cracks or chips and inspect membranes for holes, tears or warping. If you notice any damage to these parts, replace them immediately.


i Info

- Carry out all the steps with care and assemble the pump kit correctly for optimum performance.


2.1 Assembly of pump kit

Wash hands before touching breastpump, kit and breasts and avoid touching the inside of containers or lids.


Snap membrane onto valve until membranes lies completely flat on valve.


Push yellow valve securely onto inner ring of pump body.

Screw bottle onto pump body.

Insert PersonalFit™ breastshield into top opening on pump body.

Insert hard yellow adaptor end of tube into pump body.

Continued next page.

! Caution


- Wash hands thoroughly with soap and water before touching the breast and pump parts.

Questions? Visit www.medela.com or call us at 1-800-435-8316.

2. Assembly (cont.)

Warning

- Use only the power adaptor supplied with your Swing Breastpump.
- Always inspect the power adaptor prior to use for damaged, kinked, twisted, frayed or exposed wires. If damaged, immediately discontinue use and call Medela Customer Service at 1-800-435-8316.


Twists, kinks or bends


Note

- See **Section 2.3.2** Battery Compartment for instructions on the use of the Battery Compartment.
- See **Section 2.3.3** Traveling Outside the U.S. for information on using your device internationally.


2.1 Assembly of pump kit (cont.)

Insert the triangle end of tube into triangle port on motor unit.

Plug the metal end of the power cord into the power adaptor port.

Plug power adaptor into a standard electrical outlet.

2. Assembly (cont.)

Caution

- While some discomfort may be felt when first using a breastpump, using a breastpump should not cause pain. If you are unsure about breastshield sizing, please contact a healthcare professional or breastfeeding specialist who can help you get a proper fit.

2.2 Correctly fitting breastshields

Properly fitting breastshields will ensure comfort and efficient milk expression.

To determine if you have the correct size, look at the following diagram:

- The nipple should be centered in the tunnel of the shield.
- The nipple should move freely without rubbing while pumping.
- Minimal or no part of your areola should be pulled into the tunnel of the breastshield.
- You should see a gentle, rhythmic motion in the breast with each cycle of the pump.
- After pumping, your breast should feel much less full with no areas of hardness.

Correct fit

Medela offers smaller and larger sizes which are available at many local retailers or visit www.shopmedela.com.

PersonalFit Breastshields

21 mm (Small)

Item 87072

24 mm (Medium)

Item 87073 (included)

27 mm (Large)

Item 87274

30 mm (X-Large)

Item 87075

36 mm (XX-Large)

Item 87084

Choosing the right size PersonalFit™ breastshield:

- Determine the size you're currently using. If you're not sure, look for the size embossed on your breastshield (see picture). 24 mm (M) breastshields are provided with Medela breastpumps.

- While pumping, compare your fit to the images below as a sizing guideline.

Correct fit

Your breastshield is **too small**; try a larger size


Your breastshield is **too large**; try a smaller size

2.3 Powering your pump

2.3.1 Power Adaptor Care


To care for your adaptor be sure to follow these storage directions.


Do not wrap cord of power adaptor around the plug body.

Correct!

Do not unplug power adaptor by pulling on the cord.

Correct!

2. Assembly (cont.)


! Note

- Remove the batteries if they are not going to be used for an extended period of time.
- If breastpump does not turn on, make sure the batteries have been inserted properly. If this does not solve the problem refer to the troubleshooting section of this instruction booklet or contact Medela Customer Service at 1-800-435-8316.
- Each set of batteries will provide approximately 2 hours of pumping time.

2.3 Powering Options (cont.)

2.3.2 Battery Compartment

1


Use 4 AA or alkaline batteries.

4

Replace battery cover.

2

Replace all batteries at the same time.

Remove battery cover by pushing the button toward the clip.

3

Place the **NEGATIVE (-)** end of the battery against the front of the spring by pushing the battery in and down in one motion until the battery snaps into place. The **POSITIVE (+)** end of the battery should be near the positive sign on the battery pack.

2.3.3 Traveling outside of the U.S. or Canada

While travelling internationally, we recommend purchasing a Universal Power Plug Adapter (not included). Please consult with the country you are visiting to find out what type of adaptor will work best.

For information about traveling with your breastpump and expressed breastmilk, visit www.TSA.gov.

3. Pumping

Caution

- Always inspect breastshield, pump body, white membrane, yellow valve, and tube prior to use for cleanliness. If dirty, see **Section 1.2**.
- If you experience pain while pumping seek advice from your healthcare professional or breastfeeding specialist. See **Section 6** for more information.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Make sure tube is not kinked or pinched while pumping.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful.


Info

- Research has shown that you will pump more efficiently – get more milk in less time – when pumping at Maximum Comfort Vacuum during the Expression Phase.
- You should reassess your Maximum Comfort Vacuum throughout your pumping experience because it can change during the different stages of lactation.
- Do not tilt or overfill containers when pumping.
- Do not kink tubing.

3.1 Device Operation


Wash hands before touching breastpump, kit and breasts and avoid touching the inside of containers or lids.

Switch on the breastpump with . Your pumping session will start in the Stimulation Phase.

Overview of LEDs

lights up	On
does not light up	Off
even, regular flashing	Stimulation phase
light constantly on	Expression phase
rapid flashing	Min./Max. vacuum


Center the assembled breastshield over your nipple.


After two minutes, the pump will automatically change to the Expression Phase.

- If you let-down (your milk starts to flow) sooner than 2 minutes, press the let-down button.
- If let-down has not happened during the Stimulation Phase and the pump has switched to the Expression Phase, press the let-down button once to return to the Stimulation Phase. Your pumping session will start in the Stimulation Phase.

Hold the breastshield onto your breast with your thumb and index finger. Support your breast with the palm of your hand.

The vacuum can be adjusted in every phase. Comfort vacuum: Use to increase the vacuum until it is slightly uncomfortable. Then press to dial it back minimally.

Note: The LED light will flash rapidly when you have reached minimum and maximum pumping settings.

Continued next page.

3. Pumping (cont.)

! Note

- Always disconnect Swing™ from the power source after expression.
- Only fill the bottle to the 150 mL mark.
- Swing switches off automatically 30 minutes after last setting adjustment.


i Info

- **How often should you pump?** A breastpump is a replacement for when you are separated from your baby. It is important to pump when the baby would be breastfeeding, which usually means a working mother pumps 2-3 times during an 8 hour working day.
- **How long should your pumping session last?** Pumping times can vary from mother to mother sometimes 15 minutes, sometimes up to 30 minutes.
- **How should your breasts feel after pumping?** After pumping, your breasts should feel soft and there should be no firm areas. Before pumping, your breasts will have a firm, heavy feeling. After pumping make sure your breasts are not firm or have lumpy areas. This could indicate that the breast is not draining all over. See **Section 6** for more information.


3.1 Device Operation (cont.)


Switch off the breastpump with .

Close the bottle with a lid. Follow the instructions in **Section 5** "Storing and thawing breastmilk". Use bottle stand to prevent the bottle from tipping over.

Unplug breastpump from power source.

Disassemble and clean per Cleaning Instructions in **Section 1**.

3. Pumping (cont.)

Swing™ is small and lightweight, and comes with everything you need.

3.2 Mobility

Use on any tabletop or flat surface.

Clip onto your belt or pocket.

Hang around your neck, over your shoulder, or from a hook.

4. Storing, Preparing and Feeding Breastmilk

4.1 Storing breastmilk

Warning

- Do not thaw frozen breastmilk in a microwave or in a pan of boiling water.
- Do not microwave breastmilk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breastmilk.)
- When freezing, do not fill container more than 75% full to allow space for expansion. Label container with date of pumping.

Freshly Expressed Breastmilk Storage Guidelines (For Healthy Term Babies)

Room Temperature	Cooler with Ice Pack	Refrigerator	Freezer	Thawed Breastmilk
4–6 hours at 66–78 °F (19–26 °C)	24 hours at 59 °F (15 °C)	3–8 days at 39 °F or lower (4 °C)	6–12 months 0–4 °F (-18–20 °C)	use within 24 hrs

References: www.BreastmilkGuidelines.com

4.2 Preparing breastmilk

- Thaw breastmilk overnight in the refrigerator. Thawed breastmilk is safe in the refrigerator for 24 hours. Do not refreeze thawed breastmilk.
- Quickly thaw breastmilk by holding the bottle under warm running water.
- Place the sealed bottle in a bowl of warm water for 20 minutes to bring it to body temperature.
- If you are adding expressed breastmilk to a container of already frozen breastmilk, make sure to add a lesser amount than the already frozen amount.

4.3 Feeding breastmilk


It is recommended that breastfeeding is well established prior to bottle feeding your baby.

- Always inspect the bottle, nipple and other components immediately before and after each use. If nipple appears cracked or torn, discontinue use immediately.
- To prevent possible choking hazard, test strength of nipple by pulling on bulb portion of the nipple.
- Do not attempt to enlarge the nipple hole.
- Infants must not be bottle fed without adult supervision.
- Nipple should not be used as a pacifier.

5. Troubleshooting

Low or No Suction	Valve and Membrane Check - Remove and separate valve from membrane. - Inspect the valve for chips or cracks. Inspect the membrane for holes or tears. If flawed or damaged replace valve and/or membrane prior to pumping. - Wash valve and membrane per Cleaning section. - Reassemble. Make sure the membrane lays flat on the valve, and if it does not, replace the valve and/or membrane prior to pumping. Breastshield and Pump Body Check - Inspect breastshield and pump body to make sure there are no cracks and they securely attach to each other. - Make sure the breastshield forms a complete seal around the breast. Connection Point Check Inspect following connection points to ensure attachment is secure - Tubing adapter to pump body. - Tubing end to motor unit. - Check that the tubing is not kinked or pinched. - Other: Perform steps of troubleshooting Section "milk in motor". - Test vacuum. If suction is not improved, contact Medela Customer Service at 1-800-435-8316.
Power source	- Make sure pump is plugged in correctly. Pump will NOT turn on if tubing is plugged into the AC adapter port. - Make sure the "on" button is being held between 1 and 3 seconds. - Do not repeatedly hit the on/off button, or hold the on/off button down longer than 3 seconds. - If using batteries, make sure transformer is not plugged into pump. - If using batteries, make sure batteries are installed correctly inside pump. - If pump is still not working with the power adaptor or the battery compartment, contact Medela Customer Service at 1-800-435-8316.
Condensation or milk in tubing	Refer to "Tubing and motor unit cleaning" section 1.4.
No "let-down" or breastmilk expression	- Ensure your breastpump is assembled correctly and that there is suction. - Relax and take a 10-15 minute break if let down is not achieved after 2 consecutive pumping sessions. - Consult with your healthcare provider or breastfeeding specialist if expression does not occur.

Milk in Motor

Turn off the pump and unplug from power source.

- Open the battery cover, remove the expanded-plastic liner. Wash the expanded-plastic liner with clear, lukewarm water and squeeze out well. Allow to dry thoroughly.
- Pour about 2 oz of warm soapy water into a container.
- Connect the tubing and the power cable onto the pump.
- Hold the device with the control button upwards over a basin and switch on. Press  and ensure that the lowest vacuum level is set by pressing  twice.
- Hold the end of the tubing that is not connected to the pump in the container with the soapy water.
- The water is pumped through the tubing into the device. Allow it to pump until all the water has been flushed through.
- Repeat procedure with 1 oz of clear, warm water.
- Then leave the device switched on for 3 minutes, until the pump is completely dry.
- As soon as the expanded-plastic liner is dry, put it in place again and close the battery cover.
- If necessary, clean the pump set according to cleaning instructions "Prior to first use and once daily".

6. Supplemental Information/Disposal

Supplemental Information

Following are some common breastfeeding related conditions. If you experience any of these symptoms, contact a healthcare professional or breastfeeding specialist.

	Symptom	Potential Causes
Engorgement	Breasts are hard and uncomfortable, possibly with reddened areas.	Milk not draining properly from breast. Milk can drain back into the tissue causing areas to swell and become tender.
Blocked/Plugged Ducts	Area of breast will look reddened and may be tender to touch.	Milk not being drained from a specific duct. The area becomes 'clogged' up and milk is then prevented from flowing.
Mastitis	It often follows on from engorgement but can just occur suddenly. The initial symptoms are similar to flu – fatigue, headache and muscle aches, fever and localized breast tenderness. It can be in one or even both breasts and needs immediate treatment.	A bacterial infection in the breast tissue which is often accompanied by cracked nipples.

A Healthcare Professional or Breastfeeding Specialist can give you guidance in breastfeeding your baby and using a breastpump. Their guidance is valuable in terms of the lifetime of health benefits for both you and your baby. A listing of Lactation Professionals are available at www.medelalocator.com or call 1-800-TELL-YOU (1-800-835-5968).

If you have not resolved the problem with your breastpump or you have further questions, please contact Medela Customer Service at 1-800-435-8316 or email customer.service@medela.com.

Disposal Guideline

Properly dispose of the device in accordance with local, state and federal regulations.

7. Warranty

Caution

This is a single user product. Use by more than one person may present a health risk and performance questions and voids the warranty.

This product is warranted by Medela, Inc. to the original retail purchaser to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and accessories) from the date of purchase. In the event of a defect, Medela will repair or, at Medela's option, replace this product, without charge for such replacement, parts or labor. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO A DURATION OF 90 DAYS FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

Before you make a claim under this warranty, it may save time and expense to call Medela Customer Service (toll free) at 1-800-435-8316. You may also call this number for additional information concerning this warranty. If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, prepaid, together with your dated bill of sale or other proof of purchase and a brief statement of the problem to the following address:

Medela, Inc. – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS

Medela Canada, Inc.,
4160 Sladeview Crescent, Unit #8,
Mississauga, Ontario,
Canada L5L 0A1
ATTENTION: RETURNS

Call first for authorization number. Returns not accepted without an authorization number.

Questions? Visit www.medela.com or call us at 1-800-435-8316.

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Medela, Inc. – Retours, Porte 4501