U by Kotex® compact tampons

2 1-800-335-6839 ubykotex.com

Scanned + Archived
By:
FindMyManual.org

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IMPORTANT Please read before using U by Kotex® compact tampons. Use the lowest absorbency for your needs.

ATTENTION: Important Information About Toxic Shock Syndrome. Sudden fever (usually 102°F or more), vomiting, diarrhea, fainting or near fainting when standing up, dizziness, or a rash that looks like sunburn can all be warning signs of a rare but serious illness, Toxic Shock Syndrome ("TSS"). TSS can rapidly progress from flu-like symptoms to an illness that can be fatal. Should these symptoms occur, discontinue use and consult a physician immediately. You should also consult a physician before using tampons if you have had TSS warning signs in the past.

Women using tampons during their menstrual period are susceptible to the risk of contracting TSS. The reported risk of TSS is higher among teenage girls and women under 30 years of age, but it can occur at any age. The incidence of TSS is estimated to be 1 to 17 per 100,000 menstruating women and girls per year.

Studies indicate that higher absorbency tampons increase the risk of contracting TSS; therefore, we suggest you use tampons with the minimum absorbency needed to control menstrual flow in order to reduce the risk of contracting TSS. You may avoid the risk of getting tampon-associated TSS by not using tampons, and reduce the risk by alternating tampon use with the use of pads or pantliners during your menstrual period. Please consult a physician if you have any further questions.

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USAGE INSTRUCTIONS Getting Ready

- After washing hands, remove tampon from wrapper by tearing along the dashed line.

Extending the Applicator

- Extend the compact applicator by pulling on the smaller inner tube while holding the textured grip area. Pull firmly until it locks into place.
- Gently tug on the string to make sure it is hanging outside the tube.

Insert the Tampon

- Grasp textured grip, keeping pointer finger free.
- Relax and find a comfortable position, either sitting/hovering over the toilet, or standing with one foot elevated.
- Insert the applicator into your vagina at a slight upward angle, until your fingers are just touching your body. With your pointer finger, push the inner tube all the way into the outer tube, with the string outside your body.
- Withdraw the applicator and dispose.
- Do not flush the applicator or wrapper. Wash your hands.
- The tampon should feel comfortable. If you feel any discomfort, the tampon may not be placed properly – remove it and insert a new one.

Removing the Tampon

- Relax and pull the string at the same angle as the tampon was inserted. Dispose in a trash receptacle. DO NOT FLUSH tampon, applicator, or wrapper.
- Change your tampon at least every 4 to 8 hours. Remove each tampon before you insert a new one. Remove the last tampon at the end of your period.
- If you have trouble getting the tampon out, see your doctor or clinic for help.

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