HOMEDICS®

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ParaSpa™ Plus

Paraffin Bath

Scanned + Archived By: FindMyManual.org

Please Take A Moment Now

Register

Your Product At:
www.homedics.com/register

Your valuable input regarding this product will help us create the products you will want in the future.

Distributed by
HOMEDICS®

HoMedics, LLC
3000 Pontiac Trail
Commerce Township, MI
48390

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Instruction Manual and Warranty Information

PAR-350
PAR-350-THP

2 year limited warranty

HoMEDICS®

IMPORTANT SAFEGUARDS

WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY WHEN CHILDREN ARE PRESENT, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:

READ ALL INSTRUCTIONS BEFORE USING

DANGER – TO REDUCE THE RISK OF ELECTRIC SHOCK:

- ALWAYS unplug this appliance from the electrical outlet immediately after using and before cleaning.
- DO NOT place or store appliance where it can fall or be pulled into a tub or sink.
- DO NOT place in or drop into water or other liquid.
- DO NOT reach for a product that has fallen into water. Unplug immediately.
- DO NOT use while bathing or in a shower.

WARNING – TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK, FIRE OR INJURY TO PERSONS:

- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
- DO NOT operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
- Use this appliance only for its intended use as described in this manual. DO NOT use attachments not recommended by HoMedics.
- NEVER operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the appliance to a HoMedics Service Center for examination and repair.
- DO NOT carry this appliance by supply cord or use cord as a handle.
- Keep the cord away from heated surfaces.
- NEVER operate on a soft surface such as a bed or couch where the air openings may be blocked.
- NEVER drop or insert any object into any opening.
- DO NOT use outdoors.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- To disconnect, turn all controls to the off position, then remove plug from outlet.
- DO NOT handle plug with wet hands.
- Use heated surfaces carefully. May cause serious burns. DO NOT use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
- NEVER use while sleepy or drowsy

This unit is designed for household use only. This unit is not designed for commercial use.

Caution: All servicing of this product must be performed by authorized HoMedics service personnel only.

SAVE THESE INSTRUCTIONS

CAUTION—PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE OPERATING.

SPECIAL WARNINGS:

IF YOU HAVE OR HAVE HAD ANY SKIN CONDITIONS, DISEASES AND OR HEALTH PROBLEMS, YOU SHOULD CONSULT YOUR PHYSICIAN BEFORE USE.

- This product is not intended for medical use.
- Not recommended for use by people with diabetes.
- Use bath only when seated. DO NOT stand on or in the tub unit.
- DO NOT apply paraffin wax on open cuts or wounds, warts, inflamed skin, or when the sensation of any extremities is decreased.
- Consult your physician if you have any questions or concerns about acute inflammation and peripheral vascular disease or decreased sensation of any extremity before using.
- Discontinue use and consult a physician if the paraffin wax feels too warm or if it does not feel warm enough. This may indicate a possible physical problem.
- Discontinue use and consult your physician if dermatitis, discomfort, irritation, or pain occurs from continued paraffin treatments.
- Test the paraffin wax on the inside of your wrist to make sure the temperature is comfortable.
- DO NOT overfill the unit and always keep the paraffin wax level about 2 inches from the top of the bath.
- NEVER operate the unit without paraffin wax.
- DO NOT melt paraffin wax in its plastic packaging, nor melt it over an open flame, on a range, stove or in a microwave.
- DO NOT pour the melted paraffin wax from the tub/unit.
- Unplug unit before filling with paraffin wax.
- DO NOT use accessories or wax not recommended or sold by HoMedics.
- If the supply cord is damaged, it must be replaced by HoMedics' Service Center.
- DO NOT put water into the tub/unit.
- DO NOT add oil or any other additives to the paraffin wax.
- This appliance is intended for household use only.
- **Attention California Residents:**
 WARNING: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

Instructions for use:

This appliance has a polarized plug (one blade is wider than the other). To reduce the risk of electric shock, the plug is intended to fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not change the plug in any way.

This unit is controlled by a thermostat. The unit can hold a maximum of 3 lbs. of paraffin wax for operation. The melting time is less than 2 hours and it is designed to run continuously for 24 hours a day.


IMPORTANT NOTE:
This unit may feel very warm to sensitive skin. However, in order to fully experience the benefits of warm paraffin, the wax must be melted to a temperature between 123–131 degrees Fahrenheit. For proper application, quickly dip the body area (i.e. hand or elbow) into the melted paraffin wax and remove immediately. Wait a few seconds for wax to solidify. This should be repeated 6–10 times. Never soak in the wax for more than 1–2 seconds.

If you are very sensitive to heat or feel that this appliance may not be for you, please return this product to the place of purchase, **prior to use.**

SERVICING OF DOUBLE-INSULATED APPLIANCES
In a double-insulated appliance, two systems of insulation are provided instead of grounding. No grounding means is provided on a double-insulated appliance, nor should a means for grounding be added to the appliance. Servicing a double-insulated appliance requires extreme care and knowledge of the system, and should be done only by qualified service personnel. Replacement parts for a double-insulated appliance must be identical to the parts they replace. A double-insulated appliance is marked with the words "DOUBLE INSULATION" or "DOUBLE INSULATED." The symbol "⧈" may also be marked on the appliance.

How to set up the ParaSpa Plus unit:

1. Remove all contents from packaging and check to make sure you have all of the following components of the ParaSpa Plus paraffin bath:
 - Paraffin bath tub with locking lid
 - Plastic grate
 - Three 1 lb. bags of Paraffin Wax
 - 20 plastic liners
2. Place unit on a level and sturdy table or surface and open the lid.
3. Make sure ParaSpa Plus tub is empty before putting wax into tub.
4. Open the bags of Paraffin Wax and place wax in the tub. Open only as many bags as you plan to use. The unit is designed to hold a maximum of 3 lbs. of melted paraffin wax and a minimum of 2.5 lbs (Fig. 1). Always keep the level of melted paraffin wax filled to the "MIN" level marked on the inside of the tub.
5. Plug the unit into a 120V A/C household outlet.
6. Close the locking lid. Push power button. "Heating" light illuminates red while wax is melting (Fig. 2). Most of the paraffin wax should melt in 3 hours.
7. "Ready" light will come on when wax is completely melted. Open the lid and carefully place the grate into the wax and allow it to sink to the bottom of the tub (Fig. 3). The ParaSpa Plus unit is designed to be used at room temperature. If the wax starts to solidify or flakes form on the top surface, close lid and press power button.
 NOTE: Initial wax melt time will take approximately 3 hours. If green ready light comes on before wax is fully melted, replace lid and wait until 3 hours has been reached for complete melt.

Fig. 1




Fig. 2

Fig. 3

NOTE:
To avoid burns, be sure to add the grate before dipping.

Fig. 4

Fig. 5

NOTE: The temperature of the room can have an effect on the unit's performance. Therefore we always recommend testing the paraffin temperature on the inside of the wrist area before dipping your entire hand or elbow.
NOTE: NEVER open lid or dip your hand, foot or elbow while the "heating" light is flashing.
NOTE: Place the grate into the unit only when the paraffin wax is melted.

How to use the paraffin bath:

1. Before immersing your hand or elbow test the paraffin wax on the inside of your wrist to make sure the temperature of the wax is comfortable. Remove all jewelry before dipping (i.e. rings, bracelets, watch).
 NOTE: It is suggested you wash the area to be treated and apply lanolin-based hand lotion over the area before dipping it into the paraffin wax. This will allow the wax to be removed easier and will help moisturize the skin.
2. Dip your hand, elbow or foot into the melted paraffin wax and remove immediately. Wait a few seconds until the wax solidifies.
 Try to avoid moving your hand or elbow while wax is applied.
 NOTE: Dipping process is quick. Do not leave hand or elbow in the tub. Do not soak. Dip quickly in and out.
3. Repeat this step several times until several layers have formed; 6–10 coatings are suggested.
4. For extended warmth, place a plastic liner (provided) over the treated area and twist bag at opening. Keep paraffin wax and liner on treated area for 10–20 minutes for ultimate luxury (Fig. 4).
5. Peel off the wax from the treated area. If using a plastic liner, peel the plastic lining from the area being treated. Paraffin wax will peel off with liner (Fig. 5).
6. Discard used wax.
 NOTE: If you are the only person using the ParaSpa Plus, the wax can be re-used. Simply put the wax back into the tub to re-melt.


Maintenance:


Periodically, the used portion of paraffin wax should be replaced and the unit should be cleaned.

1. Unplug unit and remove grate. Allow paraffin wax to solidify overnight.
2. Push on one end of paraffin wax block. This will tip it up. Lift the block out of the unit and dispose of it. DO NOT USE SHARP OBJECTS TO REMOVE THE PARAFFIN WAX BLOCK, THIS MAY SCRATCH OR DAMAGE THE TUB. If the paraffin wax does not remove easily, plug unit in and allow to heat for 10–15 minutes, then remove the paraffin wax block and unplug the unit.
3. Wipe the inside of the unit with a disposable towel to remove residual dirt and paraffin.
4. With a damp, soapy cloth, wipe clean. DO NOT POUR WATER INTO THE UNIT, OR IMMERSE THE UNIT IN WATER.
5. Store in a cool, dry place. When you are ready to re-use the ParaSpa Plus paraffin bath for the next heat treatment, repeat steps 2–7, "How to set up the ParaSpa Plus unit". Then, follow steps 1–6, "How to use the paraffin treatment."

ParaSpa™ Plus

Paraffin Bath

20 Plastic Liners

3-1lb. Bags of Paraffin Wax

HOMEDICS®

TWO YEAR LIMITED WARRANTY

HoMedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of two years from the date of original purchase, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your HoMedics product, contact a consumer relations representative for assistance. Please make sure to have the model number of the product available.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of HoMedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from country to country. Because of individual country regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com. For Canada, please visit: www.homedics.ca.

For Service in the USA:

Email:
cservice@homedics.com

Monday-Friday
8:30am-7:00pm EST

1.800.466.3342

For Service in Canada:

Email:
cservice@homedicsgroup.ca

Monday-Friday
8:30am-5:00pm EST

1.888.225.7378


HoMedics®, TheraP, and ParaSpa are trademarks of HoMedics, LLC
©2012–2015 HoMedics, LLC. All Rights Reserved.

IB-PAR350D

Distribuido por
HOMEDICS®

ParaSpa™ Plus

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Distribuido por
HOMEDICS®

HoMedics, LLC
3000 Pontiac Trail
Commerce Township, MI
48390

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